PCI Guide Wire SION series

SYMBOLS	1
English [EN]	4
Dansk [DA]	7
Nederlands [NL]	10
Suomi [FI]	14
Français [FR]	17
Deutsch [DE]	20
Ελληνικά [EL]	24
Magyar [HU]	28
Italiano [IT]	31
Polski [PL]	35
Português [PT]	39
Español [ES]	42
Svenska [SV]	46
Česky [CS]	49
Norsk [NO]	52
Slovenský [SK]	55
Eesti keel [ET]	58
Latviešu [LV]	61
Lietuvių kalba [LT]	64
Türkçe [TR]	67
Русский [RU]	70
Български [BG]	74
Hrvatski [HR]	78
Română [RO]	81
Slovenščina [SL]	85
Malti [MT]	88

SYMBOLS

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PCI Guide Wire SION series

INSTRUCTIONS FOR USE

Read these instructions carefully before using this guide wire and observe the Indications for Use, Contraindications, Warnings, Precautions, and How to Use sections in this Instructions for Use. Failure to do so may result in complications, including serious injury to the patient or death.

If a user and/or a patient comes to know a serious incident that has occurred in relation to the guide wire, report the fact to the manufacturer and the competent authority of the Member State.

These Instructions for Use apply to the PCI Guide Wire SION series. For details (length of the guide wire, length of radiopaque section, etc.), refer to the product label.

Descriptions

PCI Guide Wire SION series has been developed as a soft-tip guide wire to reduce the risk of vascular perforation when held in a vessel, as holding a guide wire with a hard tip in a vessel may result in vascular perforation.

The basic structure of SION consists of a taper core wire, two types of distal coils (an outer coil and an inner coil) and a safety wire. The taper core wire, the coils and the safety wire are soldered. A partial length of the outer coil is radiopaque so that its position can be easily confirmed under radioscapy. SION has a flexible tip load of 0.7gf* to ensure safety in the blood vessel. An inner coil made of rope is placed inside the outer coil, and the safety wire is bonded to the inner coil and the tip so that the safety wire will not come off, even in the unlikely event that the taper core wire comes off the tip, preventing the outer coil and taper core wire from coming off. The distal side of SION is coated with hydrophilic polymer.

SION blue has the same basic structure as SION, but some of the coating on the outer coil tip is silicone coating, making the tip more flexible than that of SION, and the tip load is 0.5gf*.

The basic structure of the distal side of SION blue ES is the same as that of SION blue and the tip load is 0.5gf*, but the outer coil length is shorter and the shaft stiffness is higher. The internal structure does not have a safety wire inside, but the distal design of the taper core wire is devised so that the proximal side of the guide wire is not easily dislodged. As with SION blue, the coating of SION blue ES is silicone coating at the distal and hydrophilic at the proximal, and the hydrophilic coating length is longer than that of SION blue.

SION black differs from SION and SION blue in that it has a polyurethane resin coating on the outer coil, and a hydrophilic coating on top of that. The tip load is 0.8gf*.

The basic structure of SUOH 03 is the same as those of SION blue and SION blue ES, but it uses an outer coil made of rope, resulting in the lowest tip load in the SION series at 0.3gf*. The distal surface is coated with hydrophilic coating.

In all PCI Guide Wire SION series, the proximal side of the taper core wire is coated with PTFE (polytetrafluoroethylene). The nominal outer diameter of the product is 0.36 mm (0.014 inch) over the complete length. At the proximal end of the devices less than 300 cm, there is a connection for an extension wire, which is used to replace the combined device; ASAHI INTECC's Extension Wire can be connected.

	Tip Load*	Outer Diameter	Coating Length	Coil Length	Radiopaque Length
SION	0.7gf	0.36mm (0.014")	Hydrophilic coating : 28.5cm	28cm	3cm
SION blue	0.5gf	0.36mm (0.014")	Silicone coating : 1.5cm Hydrophilic coating : 19.0cm	20cm	3cm
SION blue ES	0.5gf	0.36mm (0.014")	Silicone coating : 2.5cm Hydrophilic coating : 39.5cm	8.5cm	3cm
SION black	0.8gf	0.36mm (0.014")	Polyurethane coating : 20.0cm Hydrophilic coating : 40.0cm	12cm	3cm
SUOH 03	0.3gf	0.36mm (0.014")	Hydrophilic coating : 52.0cm	19cm	3cm

Indications for Use

This guide wire is intended for general percutaneous coronary intervention (PCI). This guide wire is intended to facilitate the placement of balloon dilatation catheters or other interventional therapeutic devices within cardiovascular. The average duration of use is expected to be 120 minutes.

This guide wire is not intended for use in the neurovasculature and peripheral vasculature.

Target Treated Population

PCI Guide Wire SION series are intended for patients who meet indications for PCI procedures. It is not indicated for use on neonates, infants, or children as well as patients who are pregnant, have impaired renal function, and who have an allergic reaction to contrast media or polyethylene glycol (PEG).

Contraindications

- Confirm the compatibility of guide wire diameter with other interventional device before actual use.
- Do not use the guide wire in patients who have an allergic reaction to contrast media or who have impaired renal function.
- This guide wire contains PEG400 in its hydrophilic coating. Do not use the guide wire in patients who have allergic reactions to polyethylene glycol (PEG).
- Do not use the guide wire for pregnant patients or patients who may be pregnant (X-ray may affect the fetus).
- This guide wire must be used in an institution where emergency surgical operation can be performed immediately. If an emergency surgical operation is unavailable, in the worst case, life-threatening events may occur.
- This guide wire is presterilized with ethylene oxide gas (EOG) and is intended for single use only. Do not reuse or resterilize. If reused or resterilized, the performance or quality of this guide wire and package may be compromised and there is a risk of complications, including infection.

Warnings

- This guide wire must be used only by a physician who is fully trained in PCI treatment.
- Never use this guide wire in combination with catheters having metal parts that may directly contact the surface of this guide wire (e.g. atherectomy catheters, metal dilators), metal needles or metal sheath inserters.
- During use of the guide wire, check and monitor the movement of its tip under fluoroscopy. Do

not use in areas of vessel that are not or cannot be visualized.

- Do not connect this guide wire with detachable extension wires produced by manufacturers other than ASAHI INTECC. Do not manipulate this

guide wire while the detachable extension wire is connected. Otherwise, the guide wire may be damaged, or the detachable extension wire may be unintentionally detached.

- Do not modify this guide wire for any reason.

Precautions

- Do not use the guide wire after the expiration date indicated on the label. Discard any guide wire that exceeds the expiration date.
- Guide wires are delicate instruments and should be handled carefully.
- If the package is opened or damaged, do not use the guide wire. Do not open the package until just prior to use. Use aseptic technique in handling and using the guide wire.
- Always keep the guide wire in a wet state to prevent damage to its hydrophilic coating.
- Do not wipe this guide wire using an organic solution such as alcohol.
- Take preventive measures against infection after use. Discard this guide wire and package as medical waste.
- Never push, auger, withdraw, or torque this guide wire that meets resistance. Torquing or pushing this guide wire against resistance may cause damage and/or tip separation of this guide wire or direct damage to a vessel. Failure to do so may cause damage to this guide wire. Determine the cause of resistance under fluoroscopy and take any necessary remedial action.
- Always advance or withdraw the guidewires slowly and carefully.
- When torquing this guide wire inside the blood vessel, do not torque continuously in the same direction. When torquing the guide wire, rotate it clockwise and counterclockwise alternately. Do not rotate more than two rotations (720°) in the same direction.
- Flush the guide wire with heparinized and sterilized saline or other suitable solution while removing and reinserting it to prevent air from entering the interventional device. Perform exchange of this guide wire carefully to prevent air entry and/or trauma. When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and is not against the vessel wall. Failure to do so may result in trauma. Use the radiopaque marker of the interventional device to confirm position.
- Free movement of the guide wire within the interventional device is an important feature of a

steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostasis valve with an adjustable valve if it is found to inhibit the guide wire movement.

- Do not perform stent placement using more than one guide wire or wire operation through stent strut. Otherwise, the stent may be damaged or the guide wire may break or break apart.
- Before use, confirm that this guide wire is appropriate in size and shape for the procedure and is compatible with the device to be used in combination.
- Refer to the Instruction for Use attached to the medical device (e.g. balloon dilatation catheters or other interventional devices) to be used in combination with this guide wire.
- When connecting the detachable extension wire with this guide wire, it should be securely inserted to the boundary line between non PTFE coating and PTFE coating at the proximal end of the guide wire. If any abnormality is found in the connection, stop using the guide wire immediately.
- Inspect the guide wire carefully for bends, kinks, or other damage prior to use and whenever possible during the procedure. Do not use the guide wire if it is damaged.
- When shaping the tip, be careful not to damage the coil.
- Before inserting the guide wire, check which end is the distal end and make sure to insert it from the distal end.
- Maintain a continuous saline flush between the guiding catheter and the interventional device and between the interventional device and the guide wire during the procedure. Flushing prevents contrast media formation and/or clotting on the guide wire and in the catheter lumen.
- For the holder with the distal end protection tube (refer to How to Use, Figure 1), do not remove or insert the distal end protection tube while the guide wire is housed in the holder.
- For the holder with the distal end protection tube (refer to How to Use, Figure 1), remove the distal

end protection tube before housing the guide wire in the holder again.

- This device contains the following substance defined as CMR 1A/1B and endocrine-disrupting

substance in a concentration above 0.1% weight by weight; Cobalt: CAS No. 7440-48-4.

Malfunctions and Adverse Events

Possible malfunctions and adverse events through usage of this guide wire include, but are not limited to the following:

1. Malfunctions

Failure to cross a lesion	Breakage or bending of the guide wire	Damage such as separation
Difficulty in removal	Coming off of coating	Debris from a damaged guide wire

2. Adverse events

Death	Infection	Vessel dissection
Vessel perforation	Bleeding complications	Embolism
Thrombus	Infarction	Ischemia
Arrhythmia	Vasospasm	Vascular occlusion
Blood pressure reduction	Allergic reaction	

How to Use

1. Inspection prior to use

- Before use, inspect carefully and confirm the guide wire and package are undamaged.
- Before use, confirm that the guide wire is compatible with the interventional devices to be used.

2. Preparation

- Select the most suitable guide wire for the affected area and remove the holder tube containing the guide wire from the sterile pack.
- Fill up the holder tube with heparinized and sterilized saline using a syringe to soak the entire device. Take note that heparinized and sterilized saline may spill out of the holder tube during at this time. For the holder with the distal end protection tube (Figure 1), a syringe can be connected to the distal end protection tube.
- Fill the lumen of the interventional device with the heparinized and sterilized saline before inserting the guide wire.
- Release the guide wire from the tail clip and

slide the guide wire slowly, push distal end of the guide wire through the holder tube(B).(Figure 2)


Figure 1

Figure 2

visible during interventional device manipulations.

■ Rapid exchange system

- Engage the guiding catheter.
- Insert the guide wire introducer into the Y connector of the guiding catheter.
- Carefully insert the guide wire tip into the introducer.
- Advance the guide wire through the guiding catheter under fluoroscopy until the guide wire tip is just proximal to the tip of the guiding catheter.
- Attach the torque device to the guide wire if necessary.
- Advance the guide wire under fluoroscopy to pass through the lesion. Confirm by angiography that the guide wire has passed the narrowed target lesion.

g) Observe the guide wire movement in the vessels.

h) Remove the guide wire torque device and the guide wire introducer.

i) Track the interventional device over the guide wire while preventing the guide wire from moving, and advance until the lesion is reached. Ensure that the guide wire distal tip and its location in the vessel are visible during interventional device manipulations.

4. Procedures to exchange the guide wire

■ Over-the-wire system

- Remove this guide wire slowly while monitoring the movement of this guide wire under fluoroscopy.
- Insert the next guide wire in accordance with the directions in this "How to Use" section.

e) When the distal end of the guide wire is extended 5 to 6 cm beyond the other end of the holder tube (or the distal end protection tube), if necessary, shape the tip in accordance with standard practice.

f) Gently grasp the guide wire which came out from the distal end of the holder tube (or the distal end protection tube), at the point as close to the holder tube as possible and pull the guide wire out slowly and carefully.

g) If resistance is felt when removing the guide wire from the holder tube, inject heparinized and sterilized saline into the holder tube again. Continue removing the guide wire after the resistance is not felt anymore. Be sure to inject enough heparinized and sterilized saline to the holder tube because hydrophilic coating may be damaged when removing the coated guide wire forcefully.

h) After pulling the guide wire out of the holder tube, inspect it to make sure that it is not damaged.

b) Advance this guide wire carefully until its tip is just proximal to the interventional device tip.

c) Engage the guiding catheter and insert the interventional device system (with guide wire) into the Y connector.

d) Advance the interventional device system through the guiding catheter until the tip of the guide wire system is just proximal to the tip of the guiding catheter.

e) Tighten the hemostasis valve of Y connector to create a seal around the interventional device. Ensure the guide wire movement is still permitted.

f) Check to make sure the guide wire moves smoothly.

g) Attach the torque device to the guide wire if necessary.

h) Advance this guide wire under fluoroscopy to pass through the lesion. Confirm by angiography that the guide wire has passed the narrowed target lesion.

i) Observe the guide wire movement in the vessels.

j) Advance the interventional device until the lesion is reached while preventing the guide wire from moving. Ensure that guide wire distal tip and its location in the vessel are

Storage Method

Do not keep the guide wire and package in a bent and/or heavily-loaded condition. This guide wire and package must be kept out of water. Store in a cool, dark, and dry place.

Expiration Date

The expiration date is indicated on the label of the guide wire package.

Contents

5 pieces / box

Liability Disclaimer

By no means shall "ASAHI INTECC CO., LTD. and its affiliated companies" (hereinafter referred to as the "Company") be liable for accidents, personal injuries, adverse events due to any improper use of the product(s) or any other use inconsistent with these instructions. In no event shall the Company be liable for any damages either (i) arising out of storage of the product(s) after the shipment from the Company or (ii) due to selection of patients, surgery techniques, or any other medical activities by the medical institution that uses the product(s).

"SION" and "SUOH" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

*Tip load values represent the mean values of product specification values measured according to the internal standard testing method assuming the force applied to the lesion by the guide wire tip.

3. Procedures for insertion

■ Over-the-wire system

- Insert the distal end of the guide wire carefully into the guide wire lumen of the interventional device.

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©2024 ASAHI INTECC CO., LTD.
ASAHI INTECC CO.,LTD.
Emergo Europe B.V.