ASAHI Gaia Next

PCI Guide Wire

SYMBOLS	1
English [EN]	4
Dansk [DA]	7
Nederlands [NL]	10
Suomi [FI]	13
Français [FR]	16
Deutsch [DE]	19
Ελληνικά [EL]	22
Magyar [HU]	25
Italiano [IT]	28
Polski [PL]	31
Português [PT]	34
Español [ES]	37
Svenska [SV]	40
Česky [CS]	43
Norsk [NO]	46
Slovenský [SK]	49
Eesti keel [ET]	52
Latviešu [LV]	55
Lietuvių kalba [LT]	58
Türkçe [TR]	61
Русский [RU]	64
Български [BG]	67
Hrvatski [HR]	70
Română [RO]	73
Slovenščina [SL]	76
Malti [MT]	79

SYMBOLS

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ASAHI Gaia Next PCI Guide Wire

INSTRUCTIONS FOR USE

Read these instructions carefully before using this guide wire and observe the Indications for Use, Contraindications, Warnings, Precautions, and How to Use sections in these Instructions for Use. Failure to do so may result in complications, including serious injury to the patient or death.

If a user or a patient comes to know a serious incident that has occurred in relation to the guide wire, report the fact to the manufacturer and the competent authority.


These Instructions for Use apply to the ASAHI Gaia Next. For details (length of the guide wire, length of radiopaque section, etc.), refer to the product label.

Descriptions

This guide wire is a steerable guide wire with a shapeable tip and is available in pre-shaped variants. The outside diameter of the guide wire is 0.36 mm (0.014 inch). This guide wire is compatible with existing balloon dilatation catheters and other interventional devices used in common percutaneous coronary intervention (PCI). The distal portion of the guide wire tip is radiopaque.

For lubricity, the distal portion of the device is coated with a hydrophilic polymer and the proximal portion of the guide wire is coated with polytetrafluoroethylene (PTFE).

At the proximal end of the device, which is less than 300 cm, there is a connection for an extension wire, which is used to replace the combined device; ASAHI INTECC's Extension Wire can be connected, and the total length after connection is 355-400 cm.

Indications for Use

This guide wire is intended for general percutaneous coronary intervention (PCI). It can be used to introduce and position balloon dilatation catheters and other interventional devices within cardiovascular system.

This guide wire is not intended for use in the neurovasculature and peripheral vasculature.

Contraindications

- Confirm the compatibility of guide wire diameter with other interventional device before actual use.
- Do not use the guide wire in patients who have an allergic reaction to contrast media or who have impaired renal dysfunction.
- This guide wire contains PEG400 in its hydrophilic coating. Do not use the guide wire in patients who have allergic reactions to polyethylene glycol (PEG).

- Do not use the guide wire for pregnant patients or patients who may be pregnant (X-ray may affect the fetus).
- This guide wire must be used in an institution where emergency surgical operation can be performed immediately. If an emergency surgical operation is unavailable, in the worst case, life-threatening events may occur.

- This guide wire is presterilized with ethylene oxide gas (EOG) and is intended for single use only. Do not reuse or resterilize. If reused or resterilized, the performance or quality of this guide wire may be compromised and there is a risk of complications, including infection.

Warnings

- This guide wire must be used only by a physician who is fully trained in PCI treatment.
- Never use this guide wire in combination with catheters having metal parts that may directly contact the surface of this guide wire (e.g. atherectomy catheters, metal dilators), metal needles or metal sheath inserters.
- During use of the guide wire, check and monitor the movement of its tip under fluoroscopy. Do

not use in areas of vessel that are not or cannot be visualized.

- Do not connect this guide wire with detachable extension wires produced by manufacturers other than ASAHI INTECC. Do not manipulate this guide wire while the detachable extension wire is connected. Otherwise, the guide wire may be damaged, or the detachable extension wire may be unintentionally detached.
- Do not modify this guide wire for any reason.

Precautions

- Do not use the guide wire after the expiration date indicated on the label. Discard any guide wire that exceeds the expiration date.
- Guide wires are delicate instruments and should be handled carefully.
- If the package is opened or damaged, do not use the guide wire. Do not open the package until just prior to use. Use aseptic technique in handling and using the guide wire.
- Always keep the guide wire in a wet state to prevent damage to its hydrophilic coating.
- Do not wipe this guide wire using an organic solution such as alcohol.
- Take preventive measures against infection after use. Discard this guide wire and package as medical waste.
- Never push, auger, withdraw, or torque this guide wire that meets resistance. Torquing or pushing this guide wire against resistance may cause damage and/or tip separation of this guide wire or direct damage to a vessel. Failure to do so may cause damage to this guide wire. Determine the cause of resistance under fluoroscopy and take any necessary remedial action.
- Always advance or withdraw the guidewires slowly and carefully.
- When torquing this guide wire inside the blood vessel, do not torque continuously in the same direction. When torquing the guide wire, rotate it

clockwise and counterclockwise alternately. Do not rotate more than two rotations (720°) in the same direction.

- Flush the guide wire with heparinized and sterilized saline or other suitable solution while removing and reinserting it to prevent air from entering the interventional device. Perform exchange of this guide wire carefully to prevent air entry and/or trauma. When reintroducing the guide wire, confirm that the interventional device tip is free within the vessel lumen and is not against the vessel wall. Failure to do so may result in trauma. Use the radiopaque marker of the interventional device to confirm position.
- Free movement of the guide wire within the interventional device is an important feature of a steerable guide wire system because it gives the user valuable tactile information. Test the system for any resistance prior to use. Adjust or replace the hemostasis valve with an adjustable valve if it is found to inhibit the guide wire movement.
- Do not perform stent placement using more than one guide wire or wire operation through stent strut. Otherwise, the stent may be damaged or the guide wire may break or break apart.
- Before use, confirm that this guide wire is appropriate in size and shape for the procedure and is compatible with the device to be used in combination.

- Refer to the Instruction for Use attached to the medical device (e.g. balloon dilatation catheters or other interventional devices) to be used in combination with this guide wire.
- When connecting the detachable extension wire with this guide wire, it should be securely inserted to the boundary line between non PTFE coating and PTFE coating at the proximal end of the guide wire. If any abnormality is found in the connection, stop using the guide wire immediately.
- Inspect the guide wire carefully for bends, kinks, or other damage prior to use and whenever

- possible during the procedure. Do not use the guide wire if it is damaged.
- When shaping the tip, be careful not to damage the coil.
- Before inserting the guide wire, check which end is the distal end and make sure to insert it from the distal end.
- Maintain a continuous saline flush between the guiding catheter and the interventional device and between the interventional device and the guide wire during the procedure. Flushing prevents contrast media formation and/or clotting on the guide wire and in the catheter lumen.

Malfunctions and Adverse effects

Possible malfunctions and adverse effects through usage of this guide wire include, but are not limited to the following:

1. Malfunctions

Failure to cross a lesion	Breakage or bending of the guide wire	Damage such as separation
Difficulty in removal	Coming off of coating	

2. Adverse effects

Death	Infection	Vessel dissection
Vessel perforation	Bleeding complications	Embolism
Thrombus	Infarction	Ischemia
Arrhythmia	Vasospasm	Vascular occlusion
Blood pressure reduction	Allergic reaction	Residual

How to Use

1. Inspection prior to use

- Before use, inspect carefully and confirm the guide wire and package are undamaged.
- Before use, confirm that the guide wire is compatible with the interventional devices to be used.

2. Preparation

- Select the most suitable guide wire for the affected area and remove the holder tube containing the guide wire from the sterile pack.
- Fill up the holder tube with heparinized and sterilized saline using a syringe to soak the entire device. Take note that heparinized and sterilized saline may spill out of the holder

tube during at this time.

- Fill the lumen of the interventional device with the heparinized and sterilized saline before inserting the guide wire.
- Release the guide wire from the tail clip and slide the guide wire slowly, push distal end of the guide wire through the holder tube(B).(Figure 1)


Figure 1

- When the distal end of the guide wire is extended more than 4cm beyond the other end of the holder tube, if necessary, shape the tip in accordance with standard practice.
- Gently grasp the guide wire which came out from the distal end of the holder tube, at the point as close to the holder tube as possible and pull the guide wire out slowly and carefully.
- If resistance is felt when removing the guide wire from the holder tube, inject heparinized and sterilized saline into the holder tube again. Continue removing the guide wire after the resistance is not felt anymore. Be sure to inject enough heparinized and sterilized saline to the holder tube because hydrophilic coating may be damaged when removing the coated guide wire forcefully.
- After pulling the guide wire out of the holder tube, inspect it to make sure that it is not damaged.

3. Procedures for insertion

■ Over-the-wire system

- Insert the distal end of the guide wire carefully into the guide wire lumen of the interventional device.
- Advance this guide wire carefully until its tip is just proximal to the interventional device tip.
- Engage the guiding catheter and insert the interventional device system (with guide wire) into the Y connector.
- Advance the interventional device system through the guiding catheter until the tip of the guide wire system is just proximal to the tip of the guiding catheter.
- Tighten the hemostasis valve of Yd connector to create a seal around the interventional device. Ensure the guide wire

movement is still permitted.

- Check to make sure the guide wire moves smoothly.
- Attach the torque device to the guide wire if necessary.
- Advance this guide wire under fluoroscopy to pass through the lesion. Confirm by angiography that the guide wire has passed the narrowed target lesion.
- Observe the guide wire movement in the vessels.
- Advance the interventional device until the lesion is reached while preventing the guide wire from moving. Ensure that guide wire distal tip and its location in the vessel are visible during interventional device manipulations.

■ Rapid exchange system

- Engage the guiding catheter.
- Insert the guide wire introducer into the Y connector of the guiding catheter.
- Carefully insert the guide wire tip into the introducer.
- Advance the guide wire through the guiding catheter under fluoroscopy until the guide wire tip is just proximal to the tip of the guiding catheter.
- Attach the torque device to the guide wire if necessary.
- Advance the guide wire under fluoroscopy to pass through the lesion. Confirm by angiography that the guide wire has passed the narrowed target lesion.
- Observe the guide wire movement in the vessels.
- Remove the guide wire torque device and the guide wire introducer.
- Track the interventional device over the guide wire while preventing the guide wire

from moving, and advance until the lesion is reached. Ensure that the guide wire distal tip and its location in the vessel are visible during interventional device manipulations.

- a) Remove this guide wire slowly while monitoring the movement of this guide wire under fluoroscopy.
- b) Insert the next guide wire in accordance with the directions in this “How to Use” section.

4. Procedures to exchange the guide wire

■ Over-the-wire system

Storage method

Do not keep the guide wire and package in a bent and/or heavily-loaded condition. This guide wire and package must be kept out of water. Store in a cool, dark, and dry place.

Expiration date

The expiration date is indicated on the label of the guide wire package.

Contents

5 pieces / box

Liability Disclaimer

By no means shall “ASAHI INTECC CO., LTD. and its affiliated companies” (hereinafter referred to as the “Company”) be liable for accidents, personal injuries, adverse effects due to any improper use of the product(s) or any other use inconsistent with these instructions. In no event shall the Company be liable for any damages either (i) arising out of storage of the product(s) after the shipment from the Company or (ii) due to selection of patients, surgery techniques, or any other medical activities by the medical institution that uses the product(s).

“ASAHI Gaia Next” is a trademark or a registered trademark of ASAHI INTECC CO.,LTD. in Japan and other countries.

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©2023 ASAHI INTECC CO., LTD.
"ASAHI Gaia Next" hija trademark jew trademark irregistrata ta' ASAHI INTECC CO.,LTD. fil-Ġappun u f'pajjiżi oħra.
Emergo Europe B.V.