
ALLEGRA™

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Please read all instructions carefully before using the product. Please comply with all warnings and precautions included in these instructions. Non-compliance may result in serious complications.

Implantation of the ALLEGRA Transcatheter Heart Valve may only be performed by physicians who have received appropriate training.

1. DEVICE DESCRIPTION

The ALLEGRA Transcatheter Heart Valve System (ALLEGRA THV System) is a Transcatheter Aortic Valve Implantation (TAVI) treatment and consists of the ALLEGRA Transcatheter Heart Valve (REF see chapter 1.1) and the IMPERIA Delivery System (REF see chapter 1.2) including the IMPERIA Loading System (see chapter 1.3).

1.1. ALLEGRA TRANSCATHTER HEART VALVE

The ALLEGRA Transcatheter Heart Valve (THV) is intended to be implanted within degenerative calcified aortic valves or in failing surgical aortic valve bioprostheses in a minimally invasive transcatheter implantation technique.

The valve has a tri-leaflet design and is constructed with six individual bovine pericardial parts, three parts for the skirt and three leaflets sewn to the skirt in a semi-lunar fashion to form the valve. The valve is positioned in a supra-annular position. The pericardial inner skirt creates a 12 mm sealing region against the native aortic annulus mitigating the risk of paravalvular leakage. The bovine pericardium used for the THV is treated to reduce calcification potential. The stent frame is made of laser-cut nitinol with good radiopacity and 6 additional radiopaque gold markers, positioned circumferentially to mark the height of the new valve plane.

The table below indicates the available sizes of the THV with catalogue numbers and the respective size ranges for treatment of both native aortic stenosis and degenerated surgical aortic valve (SAV) bioprosthesis (true inner diameter).

Model ALLEGRA Transcatheter Heart Valve (THV)	ALLEGRA 23	ALLEGRA 27	ALLEGRA 31
Catalogue / Reference number	THV-AO0023G1RE	THV-AO0027G1RE	THV-AO0031G1RE
Aortic annulus diameter [mm]* / Inner diameter aortic valve bioprosthesis [mm]*	16.5 - 21	21 - 24	24 - 27

*Perimeter derived

1.2. IMPERIA DELIVERY SYSTEM

The IMPERIA Delivery System (DS) is intended for transfemoral retrograde implantation of the ALLEGRA Transcatheter Heart Valve (THV).

The DS is compatible with a 0.035 inch (0.889 mm) guidewire and has to be inserted to the vessel through a separate 18 Fr introducer sheath. The catheter shaft has a distal short cartridge of 18 Fr (6.0 mm), which accommodates the THV. Behind the cartridge section, the catheter shaft reduces to 15 Fr (5.0 mm), which gives the necessary flexibility, stability and pushability. The usable length of the catheter shaft is 1150 mm and allows the implantation site to be reached in the target patient population. The atraumatic catheter tip and the loading anchor are radiopaque, and the tip of the catheter shaft has a radiopaque ring.

The DS uses a bottom-up deployment technique. By winding the release handle the catheter is retracted and exposes the inflow of the THV, which will start to open. After retracting the catheter up to the safety position, the THV is 80% deployed but still attached to the catheter by the three T-bars of the THV. In this configuration, the THV can be either implanted or recaptured and repositioned.

tioned or fully retrieved if required. After reaching the safety position and confirmation of a correct implant position the unlock wheel must be switched to the unlock position and the THV can be fully deployed by further retracting the catheter until it is fully opened and detached from the DS. The table below indicates the available models of the DS with catalogue numbers and the compatibility in regard of the THV size.

Model IMPERIA Delivery System	IMPERIA-S	IMPERIA-M
Catalogue / Reference number	DSL-AO18G2115S	DSL-AO18G2115M
Compatible with ALLEGRA THV	ALLEGRA 23	ALLEGRA 27 ALLEGRA 31
Outer diameter	18 Fr (6.0 mm)	18 Fr (6.0 mm)

The DS requires a sufficient access vessel size for placement of an 18 Fr introducer sheath. The vascular access must be carefully assessed in terms of calcification pattern/localization, tortuosity, and compatibility with the selected sheath. Clinical studies have shown safe use in access vessels of diameter ≥ 6.0 mm and also in smaller diameters depending on the pattern of any calcific atherosclerosis present.

1.3. IMPERIA LOADING SYSTEM

The IMPERIA Loading System (LS) is intended to load the ALLEGRA Transcatheter Heart Valve into the cartridge of the IMPERIA Delivery System.

The LS consists of a loading tool, a counter funnel, and a loading pin to orientate the pericardium inside the stent. The LS is included in the packaging of the DS.

2. INDICATION

The ALLEGRA Transcatheter Heart Valve is indicated for the treatment of severe calcified aortic valve stenosis in high risk patients with elevated surgical risk or in patients with a symptomatic degeneration of an aortic valve bioprosthesis.

The IMPERIA Delivery System is indicated for transfemoral delivery of the ALLEGRA Transcatheter Heart Valve.

2.1. TARGET PATIENT GROUP

The device was studied in patients with severe calcific aortic valve stenosis or with a symptomatic degeneration of an aortic valve bioprosthesis which were considered eligible for a TAVI treatment based on the clinical judgement of the heart team. This clinical judgement included frailty and co-morbidities of the patient. The degeneration of the aortic valve bioprosthesis in the valve-in-valve patient population resulted from valve stenosis, insufficiency, or the combination of both.

The target patient group shall present:

- Symptomatic degeneration of severe calcified aortic heart valve or failing surgical aortic bioprosthetic valve showing echocardiographically mean aortic gradient ≥ 40 mmHg or peak jet velocity ≥ 4.0 m/s or Aortic Valve Area (AVA) ≤ 1.0 cm² (or ≤ 0.6 cm²/m²) or symptomatic patients with severe bioprosthetic valve insufficiency.
- Current heart valve replacement guidelines (e.g., Guideline of European Society of Cardiology) should be considered.
- Anatomic requirements according to device specifications in chapter 1.1. and 1.2.

3. CONTRAINDICATION

The ALLEGRA THV System must not be used if the implanting physician believes its implantation would be contrary to the best interest of the patient.

3.1. ABSOLUTE CONTRAINDICATIONS

The ALLEGRA THV System is contraindicated in the presence of any of the following conditions:

- Unicuspid or bicuspid aortic valve.
- Non-calcified aortic stenosis.
- Combined aortic valve disease with predominant aortic regurgitation > 3.
- Distance between aortic valve basal plane and the orifice of the lowest coronary artery < 8 mm.
- Echocardiographic evidence of intracardiac thrombus or vegetation.
- Significant aortic disease such as severe obstructive calcification or marked tortuosity or kinking which would preclude a safe advancement of the delivery system.
- Iliofemoral vessel conditions such as severe obstructive calcification, severe tortuosity or kinking that would preclude safe placement of an 18 Fr introducer sheath or make endovascular access to the aortic valve impossible.
- Severe ventricular dysfunction with Left Ventricular Ejection Fraction (LVEF) less than 20%.
- Evidence of active endocarditis or other acute infections.
- Renal failure requiring continuous renal replacement therapy.
- Known hypersensitivity to contrast media, which cannot be adequately pre-medicated or contraindication to anticoagulant or anti-platelet medication (e.g., aspirin, heparin, clopidogrel) or to Nitinol (titanium or nickel) alloy or to bovine tissue.

3.2. RELATIVE CONTRAINDICATIONS

- Evidence of an acute myocardial infarction within the past 30 days.
- Cerebral Vascular Accident (CVA) within past 6 months.
- Evidence of active peptic ulcer or upper gastrointestinal bleeding within past 90 days.
- Severe mitral regurgitation.

4. VALVE-IN-VALVE SPECIFIC CONSIDERATIONS

Patient screening and procedural planning should be done as follows for treatment of degenerated SAV procedures:

- The inner diameter of the degenerated SAV bioprosthesis must be within a range of 16.5 mm to 27 mm.
- The failure mode of the degenerated SAV bioprosthesis should be identified prior to performing the procedure:
 - In case of structural damage of the failed SAV bioprosthesis, the benefit of the procedure must be weighed against the uncertain risk of complications.
 - In failed SAV bioprosthesis with structural damage, balloon aortic valvuloplasty should only be performed with extreme caution, using undersized balloons.
 - In cases of pannus or previous endocarditis, the inner diameter and inner orifice shape should be determined using Multi Slice Computed Tomography (MSCT) to assess the landing zone of the THV. These measurements should be used to select the proper size of the THV.
 - In case of patient prosthesis mismatch, the benefit of the procedure must be weighed against the uncertain risk of complications, keeping in mind that any pre-existing patient prosthesis mismatch cannot be alleviated by valve-in-valve implantation.
 - In case of additional paravalvular leaks of the failed SAV bioprosthesis, it must be considered that valve-in-valve implantation will not solve this issue.
 - In case of low position of the coronary ostia, especially in combination with shallow sinuses, the high risk of coronary occlusion must be considered.
 - The risk of coronary occlusion should be judged as higher for SAV bioprostheses with externally mounted leaflets, including (but not limited to) the Mitroflow, Crown, Trifecta and Ionescu-Shiley.

- Should the radiopaque landmarks of the failed SAV bioprosthesis be insufficient to safely identify the correct landing zone of the THV, then the leaflet calcifications and/or a pigtail at the base of the leaflets of the degenerated SAV bioprosthesis may serve as additional references.

5. PRE AND POST BALLOON AORTIC VALVULOPLASTY CONSIDERATIONS

Balloon Aortic Valvuloplasty (BAV) for pre-dilatation of the native annulus is needed for the placement of the ALLEGRA THV. The size and model of the pre-dilatation BAV catheter should be selected such that it results in effective expansion to allow full expansion of the THV upon implantation. The use of correctly sized balloons ensures effective pre-dilatation, thereby minimizing the risk of under expansion of the THV. For native aortic stenosis, the diameter of the BAV catheter used for pre-dilatation should be at least equal to the minimum annular diameter. For the ALLEGRA THV 31, at least a 24 mm balloon is recommended.

If valve function or sealing is impaired due to excessive calcification or incomplete expansion of the THV, a post-dilatation may improve valve function and sealing. If the physician determines that post-dilatation is appropriate, consider the following factors when selecting the balloon parameters to ensure patient safety: balloon catheter model, balloon size, balloon position, inflation pressure and patient anatomy.

Maximum recommended balloon diameter for post-dilatation of ALLEGRA THV:

Model THV	ALLEGRA 23	ALLEGRA 27	ALLEGRA 31
Maximum balloon size	22 mm	24 mm	26 mm

Information on BAV catheter selection:

- Overexpansion of the narrowest section of the THV resulted in damage to the commissural attachment points in bench testing and should be avoided.
- It is important to note that the mechanical compliance properties of the selected balloon influence the dilatation behavior.
- Compliant and semi-compliant balloons are softer and will more readily conform to the profile of the THV at lower pressures but must be inflated at pressures that preserve the profile of the THV.
- Non-compliant balloons are stiffer and will achieve the nominal diameter during inflation irrespective of the underlying annulus or THV resistance and might need to be downsized.
- For additional instructions on the use of BAV catheter devices refer to the specific balloon catheter manufacturer's labeling.

Information for failed SAV bioprosthesis:

- Pre-dilatation of a stenotic SAV bioprosthesis has not been systematically evaluated and should only be performed at the discretion of the heart team and after consideration of the failure mode of the failed SAV bioprosthesis.
- In failed SAV bioprosthesis with structural damage, balloon aortic valvuloplasty should only be performed with extreme caution, using undersized balloons.

6. WARNINGS

- This procedure should only be performed where emergency aortic valve surgery can be performed promptly.
- Correct sizing of the THV is essential to prevent paravalvular leakage, migration, or annular rupture. The THV is intended for use in patients with a native aortic annulus size ranging from 16.5 mm to 27 mm or with a failing SAV bioprosthesis with a true inner diameter ranging from 16.5 mm to 27 mm.
- Correct implantation height of the THV is essential to prevent interference with the anterior mitral leaflet, coronary obstruction, and poor performance of the THV.

- To avoid occlusion of the coronary ostia, the height from the aortic valve basal plane and the orifice of the lowest coronary artery should be > 8 mm in order to have a safety margin and not to cover the ostia with the skirt of the THV. The upper limit of the skirt is marked by the 6 radiopaque markers.
- Do not resterilize or re-use the THV, the DS or the LS. Resterilization or re-use may compromise the structural integrity of the products and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not use the products if the use by date has expired.
- Accelerated deterioration due to calcific degeneration of the THV (as with any glutaraldehyde cross-linked bioprosthesis) may occur in patients with an altered calcium metabolism.
- Overall durability, especially long-term, has not been established for the THV. Careful and continuous medical follow-up is advised so that THV related complications can be diagnosed and properly managed.
- It is recommended that all THV recipients be prophylactically treated for endocarditis to minimize the possibility of prosthetic valve infection.
- THV recipients should be maintained on anticoagulant therapy, except where contraindicated, as determined by their physician.
- Do not freeze. Always store the THV in a cool, dry place between 10°C and 38°C. Any THV that has been frozen must not be used for human implantation.
- Prior to the insertion of the DS into the introducer sheath, check if the DS is compatible with selected 18 Fr introducer sheath. The 18 Fr Performer™ introducer sheath from Cook Medical must not be used because it may cause damage to the tip of the cartridge.
- The safety and effectiveness of the THV implanted within a failed pre-existing transcatheter heart valve have not been demonstrated.
- Do not perform an implantation without an introducer sheath or a guidewire.
- Physicians without previous experience using the device shall be trained according to the manufacturer Physician Training Plan.

7. PRECAUTIONS

7.1. PRECAUTIONS PRIOR TO USE

- Before, during and post-implant handling of the ALLEGRA THV System shall be in accordance with the provided instructions.
- Loading of the THV may only be performed by physicians or hospital personnel who have received the appropriate training.
- Implantation of the THV may only be performed by physicians who have received the appropriate training.
- Do not use the THV if the warm temperature indicator is black, or the cold temperature indicator's bulb is purple.
- Use the DS only if the ethylene oxide indicator on the inner sterile bag is green.
- Do not use the DS for any other THV model than indicated in the instructions for use.
- Do not use the product if the integrity of the sterile packaging has been compromised.
- In the event of any non-matching information, do not use the product.
- Do not use the THV if the product jar is damaged or if the seal is broken or the outside of the product jar is wet.
- Do not use the THV if the level of the liquid inside the product jar is below the 120 ml mark.
- Do not use the THV if there are any signs of damage to the tissue or frame.
- Do not use the DS and loading components if there are any signs of damage.
- Do not use the DS if the sterile bags are damaged.

7.2. PRECAUTIONS DURING USE

- The outside of the product jar is not sterile and must not enter the sterile field.
- Follow aseptic techniques during device preparation and implantation.
- During the loading procedure all steps must be performed inside the loading bowl with cold sterile 0.9% saline solution (0°C to 8°C, with sterile ice cubes).
- Any THV that is dropped, damaged, or mishandled in any way must not be used for implantation.
- Do not handle or manipulate the THV with sharp or pointed objects.
- Do not touch or squeeze the THV during rinsing. No other objects should be placed in the rinsing bowls.
- Make sure not to rinse the THV in an already used rinsing solution.
- Do not load the THV without the protection wire.
- Do not use the DS if flushing of the system is not possible.
- Avoid repeated crimping cycles during the loading procedure.
- If you feel any resistance during loading, the loading procedure must be stopped, and the cause of the resistance has to be evaluated.
- Do not use any tool other than the loading pin, in order to prevent damage to the pericardium.
- Make sure that from the second flushing step the tip of the catheter is inside the loading bowl and that the tip is covered by the cold sterile 0.9% saline solution (0° to 8°C, with sterile ice cubes) always until the catheter is entirely closed.
- Make sure that the suture or particles from the heart valve holder do not remain on the THV.
- Make sure that there is no gap and no overlapping between the catheter and the tip.
- Keep the THV hydrated and de-aired by flushing the port of the catheter shaft until ready for implantation.
- In case of an implantation in a failing SAV bioprosthesis, carefully check the patient's anatomy regarding low coronary ostium location, especially in combination with shallow sinuses.
- In case of an implantation in a failing SAV bioprosthesis, ensure that there is no partially detached valve leaflet, which can occlude a coronary ostium in the aortic position.
- In case of an implantation in a failing SAV bioprosthesis, the risk of coronary occlusion is higher for SAV bioprostheses with external leaflets e.g., Mitroflow, Crown, Trifecta and Ionescu-Shiley.
- In case of an implantation in a failing SAV bioprosthesis, ensure that the guidewire crosses the central lumen of the failing SAV bioprosthesis and not the paravalvular space or through a defect in a leaflet or sewing ring of the failing SAV bioprosthesis.
- Ensure the diameter of the access vessels is adequate.
- Do not kink the catheter shaft during the loading procedure.
- Balloon aortic valvuloplasty of the native aortic valve is required prior to DS insertion to facilitate placement of the THV.
- For native aortic stenosis, the diameter of the aortic valvuloplasty balloon used for pre-dilatation should be at least equal to the minimum annular diameter. For the ALLEGRA 31, at least a 24 mm balloon is recommended.
- Balloon aortic valvuloplasty should only be performed when the THV is loaded into the delivery system and ready for implantation.
- In case of failing SAV bioprosthesis with structural damage, balloon aortic valvuloplasty should be performed with extreme caution.
- If you feel any resistance during introducing the DS, the procedure must be stopped, and the cause of the resistance has to be evaluated.
- Make sure that the catheter at the integrated sheath is stable to avoid any position change.
- Once the bottom of the THV is implanted in the aortic annulus, do not pull, or push on the delivery system to adjust the implant depth.
- After the safety position is reached, do not turn the release handle further.

- After full THV deployment, check if the T-bars are free from the catheter.
- Always perform THV recapture and repositioning procedure only under fluoroscopic guidance.
- Recapture and deployment of the THV can be attempted up to three times. If the THV has been fully recaptured three times, it must be removed from the patient and discarded.
- Stop closing the catheter if a resistance occurs. You may check and adjust the angulation between the THV and the cartridge.
- A higher implant position may increase the risk for paravalvular leakage. A lower implant position may contribute to a higher risk of conduction disorders.
- The catheter closure is recommended to be done in a straight portion of the descending aorta.
- Do not over-close the catheter as this may cause overlapping and make it difficult to remove the delivery system through the introducer sheath.
- If it is necessary to post-dilate the THV in the case of unacceptable paravalvular leakage or under-expansion, then proceed with caution as post-dilatation may damage the device integrity or cause migration of the THV.
- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure or inhalation of the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water. In the event of contact with eyes, seek immediate medical attention.

8. RESIDUAL RISKS AND UNDESIRABLE SIDE EFFECTS

Residual risks and potential undesirable side effects that may be associated with the use of the ALLEGRA THV System include, but may not be limited to:

- Acute coronary occlusion
- Acute kidney injury
- Acute myocardial infarction
- Acute renal failure
- Allergic reactions/intolerances (e.g., to contrast media)
- Aortic root injury (e.g., dissection, perforation)
- Arrhythmias including ventricular tachycardia or fibrillation extending to cardiac arrest
- Atrioventricular conduction disorders (e.g., AV-block, LBBB) which may require a permanent pacemaker
- Bleeding (hemorrhage)
- Cardiac tamponade
- Cardiogenic shock
- Cerebrovascular events such as TIA, Stroke
- Death
- Delivery system malfunction requiring the exchange of the device
- Device embolization / migration
- Emergent cardiac surgery (e.g., coronary artery bypass or heart valve replacement) or additional transcatheter intervention (e.g., second valve implantation or change to another device)
- Endocarditis
- Exacerbation of heart failure
- Hemolysis
- Hemorrhage requiring transfusion
- Hypertension or hypotension
- Infection
- Mitral valve injury
- Non-structural prosthetic valve dysfunction: paravalvular or/and central regurgitation, stenosis caused by mis-sizing and/or patient prosthesis mismatch
- Prosthetic valve thrombosis

- Sepsis
- Structural prosthetic valve damage (e. g. cusp tear, suture line disruption, stent fracture, calcification)
- Thromboembolism
- Vascular injury (e.g., dissection, perforation)

9. DIRECTIONS FOR USE

CAUTION: Before, during and post-implant handling of the ALLEGRA THV System shall be in accordance with the provided instructions.

9.1. PHYSICIAN TRAINING

The implanting physician should be experienced in port access, balloon aortic valvuloplasty, catheterization techniques, the TAVI implantation procedure, and trained on the use of the ALLEGRA THV and the IMPERIA DS.

9.2. REQUIRED EQUIPMENT

Standard cardiac catheterization lab equipment is required.

9.3. REQUIRED MATERIAL

The materials listed are delivered with the IMPERIA DS:

- 3 rinsing bowls incorporated in the tray (capacity approx. 500 ml / bowl)
- 1 loading bowl incorporated in the tray (capacity approx. 1000 ml / bowl)

The materials listed below are not delivered with the ALLEGRA THV or IMPERIA DS:

- 1 small bowl for heparinized saline flushing solution (capacity approx. 100 ml)
- 3 liters sterile 0.9% saline solution
- Approx. 50 ml heparinized sterile 0.9% saline solution
- Frozen sterile 0.9% saline solution (sterile ice cubes)
- Pre-shaped TAVI guidewire, 0.035 inch (0.889 mm), e.g., Safari², 300 cm length from Boston Scientific
- Percutaneous transluminal balloon aortic valvuloplasty catheter in various sizes
- Compatible 18 Fr introducer sheath
- Pigtail catheter
- 10 ml straight tip syringe without Luer-lock (for catheter flushing)
- Standard procedure supplies

CAUTION: Prior to the insertion of the DS into the introducer sheath, check if the DS is compatible with the selected 18 Fr introducer sheath. The 18 Fr Performer™ introducer sheath from Cook Medical must not be used because it may cause damage to the tip of the cartridge.

9.4. INSPECTION PRIOR TO USE

NOTE: Before using the products, check the temperature indicators, the ethylene oxide indicator, the validity of the serial or lot number, and the correct size of the implant and model of the DS by using the four-eye principle.

Carefully inspect the package before opening for any signs of abnormalities. Remove the device from the protective package and visually check that it is free of damages.

CAUTION: In the event of any non-matching information, do not use the product.

CAUTION: Do not use the product if the integrity of the sterile packaging has been compromised.

CAUTION: Do not use the THV if the warm temperature indicator is black, or the cold temperature indicator's bulb is purple.

CAUTION: Do not use the THV if the product jar is damaged or if the seal is broken or the outside of the product jar is wet.

CAUTION: Do not use the THV if the level of the liquid inside the product jar is below the 120 ml mark.

CAUTION: Do not use the THV if there are any signs of damage to the tissue or frame.

CAUTION: Do not use the DS and loading components if there are any signs of damage.

CAUTION: Do not use the DS if the sterile bags are damaged.

CAUTION: Use the DS only if the ethylene oxide indicator on the inner sterile bag is green.

CAUTION: Do not use the DS for any other THV model than indicated in the instructions for use.

9.5. DELIVERY SYSTEM HANDLING AND PREPARATION

Step	Procedure
1	Remove the double sterile bag from the carton box. Verify that all carton box and sterile bag labels match with respect to model, size, and lot number. CAUTION: In the event of any non-matching information, do not use the product. NOTE: Before opening the carton box, make sure that the DS model corresponds to the correct THV model.
2	Carefully inspect the entire sterile bags for damages. CAUTION: Do not use the DS if the sterile bags are damaged. CAUTION: Use the DS only if the ethylene oxide indicator on the inner sterile bag is green.
3	Open the outer sterile bag and ensure that only the inner sterile bag is introduced into the sterile field.
4	Remove all components of the DS from the inner sterile bag. This includes the tray system where the DS, loading tool, loading pin and counter funnel are clipped in. CAUTION: Do not use the DS and loading components if there are any signs of damage.
5	Prepare the DS and the product tray for the loading procedure as follows: Remove the rinsing bowls out of the tip part of the tray. Open the fixation strap between the handle and the tip tray and disconnect the handle tray and the tip tray from each other. Position the three parts of the tray as follows: handle part (contains the handle) on top, rinsing bowls in between and the tip part (contains the catheter) at the bottom. Attach the three parts of the tray using the connecting points.
6	Prepare the 3 rinsing bowls with at least 500 ml of sterile 0.9% saline solution at ambient temperature (15°C to 25°C) to thoroughly rinse the glutaraldehyde sterilant from the THV.
7	Prepare the loading bowl containing at least 1000 ml of cold, sterile 0.9% saline solution (0°C to 8°C, with sterile ice cubes) for the loading procedure.

9.6. THV HANDLING AND PREPARATION

CAUTION: Follow aseptic techniques during device preparation and implantation.

Step	Procedure
1	Remove the product jar from the carton box. Verify that all carton box and product jar labels match with respect to THV model, size, and serial number. CAUTION: In the event of any non-matching information, do not use the product. NOTE: Before opening the carton box, make sure that the THV model corresponds to the correct DS model.

Step	Procedure
2	Carefully inspect the entire product jar and seal for damage. CAUTION: Do not use the THV if the product jar is damaged or if the seal is broken or the outside of the product jar is wet. CAUTION: Do not use the THV if the level of the liquid inside the product jar is below the 120 ml mark.
3	Break the seal from the product jar and unscrew the screw cap. The content of the product jar is sterile and must be handled aseptically to prevent contamination. CAUTION: The outside of the product jar is not sterile and must not enter the sterile field.
4	Open and hold the product jar, while someone inside the sterile field takes out the transcatheter heart valve holder from the product jar. The THV must be inspected for any signs of damage to the frame or tissue. CAUTION: Do not use the THV if there are any signs of damage to the tissue or frame. CAUTION: Any THV that is dropped, damaged, or mishandled in any way must not be used for implantation.

9.7. THV RINSING PROCEDURE

NOTE: The whole rinsing and loading procedure requires approx. 10 minutes; 6 minutes for rinsing (3 x 2 minutes) and approx. 4 minutes for loading.

Step	Procedure
1	The THV needs to be rinsed to be free from residues of glutaraldehyde. Place the THV in the first bowl of sterile 0.9% saline solution. Be sure the solution completely covers the THV and holder. CAUTION: Do not handle or manipulate the THV with sharp or pointed objects.
2	With the THV and transcatheter heart valve holder submerged, gently agitate the THV back and forth for a minimum of 2 minutes by grasping the transcatheter heart valve holder. CAUTION: Do not touch or squeeze the THV during rinsing. No other objects should be placed in the rinsing bowls.
3	Repeat the 2-minute rinse in the second and third bowl for a total rinse time of at least 6 minutes. NOTE: Make sure not to rinse the THV in an already used rinsing solution.
4	Leave the THV in the final rinsing solution until needed to prevent the tissue from drying. CAUTION: The THV should be kept hydrated throughout the rest of the preparation procedure to prevent the tissue from drying.

9.8. THV LOADING PROCEDURE

CAUTION: Do not load the THV without a protection wire.

CAUTION: Do not kink the catheter shaft during the loading procedure.

CAUTION: Do not use the DS if flushing of the system is not possible.

CAUTION: If you feel any resistance during loading, the loading procedure must be stopped, and the cause of the resistance has to be evaluated.

CAUTION: Avoid repeated crimping cycles during the loading procedure.

CAUTION: During the loading procedure all steps must be performed inside the loading bowl with cold sterile 0.9% saline solution (0°C to 8°C, with sterile ice cubes).

Step	Procedure
1	Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution, and flush the catheter flush port.
2	Open the catheter by rotating the release handle until you reach the blue marker on the threaded rod. The rotation will also stop due to a resistance.
3	Rotate the unlock wheel from the lock position to the unlock position.
4	Continue to open the catheter until the release handle touches the rear handle. The loading anchor is fully visible, and the DS is prepared to attach the THV.
5	Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution, and flush the catheter flush port. NOTE: Make sure that during flushing the tip is in a higher position than the handle.
6	Remove the loading plate from the guiding pins. Place the THV on the loading tool, with the T-bars pointing towards inside the loading funnel. Pass the guiding pins of the loading tool through the loading guides of the heart valve holder. Mount the loading plate on the guiding pins. NOTE: For easier handling, keep the jar fixation of the transcatheter heart valve holder parallel to the THV.
7	Introduce the tip of the catheter into the guiding funnel of the loading tool. Press the THV with the loading plate against the loading funnel to carefully align the three T-bars to the grooves of the loading anchor and guide them into place. The T-bars must stay in the grooves. NOTE: Make sure not to manipulate the leaflets of the THV.
8	Rotate the release handle forward until the safety position (blue marker on the threaded rod) is reached. The three T-bars of the stent are secured by the cartridge. Continue the movement of the catheter shaft forward by turning the release handle until the commissural posts of the THV are covered by the catheter.
9	Remove the loading plate and move the remaining loading tool into the direction of the handle. It is no longer needed and must be removed after the loading process is finished. Cut the three sutures of the heart valve holder and remove the holder. CAUTION: Make sure that the suture or particles from the heart valve holder do not remain on the THV. NOTE: Take care not to manipulate the THV whilst cutting the sutures.
10	Cool the THV inside the loading bowl with cold sterile 0.9% saline solution (0° to 8°C, with sterile ice cubes). Pre-crimp the THV into the conical part of the counter funnel. Remove the counter funnel after pre-crimping.
11	Using the loading pin, push the pericardium, which protrudes outside the stent rhombs, gently back inside the stent frame. CAUTION: Do not use any tool other than the loading pin, in order to prevent damage to the pericardium. NOTE: Both ends of the loading pin can be used for the repositioning of the pericardium.

Step	Procedure
12	Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution, and flush the catheter flush port in two steps. For the first flushing step, make sure, that the tip is in a higher position than the handle. In a second flushing step, flush with the tip inside the loading bowl. CAUTION: Make sure that from the second flushing step the tip of the catheter is inside the loading bowl and that the tip is covered by the cold sterile 0.9% saline solution (0° to 8°C, with sterile ice cubes) always until the catheter is entirely closed.
13	Cool the THV inside the loading bowl with cold sterile 0.9% saline solution (0° to 8°C, with sterile ice cubes). Crimp the THV completely into the cylindrical part of the counter funnel. Do not remove the counter funnel.
14	Move the catheter shaft forward by turning the release handle until the catheter is completely closed. Close the catheter inside the loading bowl with cold sterile 0.9% saline solution (0°C to 8°C, sterile ice cubes needed). Remove the counter funnel. CAUTION: Make sure that there is no gap and no overlapping between the catheter and the tip. NOTE: If there are any air bubbles visible inside the counter funnel during the closing of the catheter, remove the counter funnel to get rid of the air bubbles.
15	Take the protection wire out of the DS.
16	Use a 10 ml syringe, filled with heparinized sterile 0.9% saline solution, and flush the guidewire lumen and the integrated sheath.
17	Activate the hydrophilic coating of the catheter carefully by wetting it with 0.9% saline solution and a sterile compress.
18	Keep the loaded DS with catheter closed in a cold 0.9% saline solution (0°C to 8°C, with sterile ice cubes) until implantation. CAUTION: Keep the THV hydrated and de-aired by flushing the port of the catheter shaft until ready for implantation.

9.9. THV DELIVERY

NOTE: Use and monitor systemic anticoagulation during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, consider alternative anticoagulation therapy.

In the case of an implantation in a failing SAV bioprosthesis:

CAUTION: Carefully check the patient's anatomy regarding low coronary ostium location, especially in combination with shallow sinuses.

CAUTION: Ensure that there is no partially detached valve leaflet, which can occlude a coronary ostium in the aortic position.

CAUTION: The risk of coronary occlusion is higher for SAV bioprostheses with external leaflets e.g., Mitroflow, Crown, Trifecta and Ionescu-Shiley.

CAUTION: During the guidewire crosses the central lumen of the failing SAV bioprosthesis and not the paravalvular space or through a defect in a leaflet or sewing ring of the failing SAV bioprosthesis.

Step	Procedure
1	The patient requires a contra and ipsilateral port with an introducer sheath in one side and with the DS compatible sheath at the other. Prepare the vascular access site according to standard technique. CAUTION: Ensure the diameter of the access vessels is adequate.
2	Introduce a pigtail catheter under fluoroscopic guidance through the introducer sheath and position it in the right or non-coronary cusp.
3	Perform aortic root angiography and find the best view of the aortic valve to acquire a position where the three leaflets (noncoronary cusp, left coronary cusp, right coronary cusp) are aligned.
4	Cross the aortic valve into the left ventricle using an angiographic catheter in combination with a 0.035 inch (0.889 mm) guidewire with a soft tip (straight or J), ideally between the non-coronary and the right coronary cusps.
5	Exchange the guidewire for a pre-shaped TAVI guidewire.
6	Remove the catheter, leaving the guidewire in place in the left ventricle, and advance the balloon aortic valvuloplasty catheter over the guidewire, cross the aortic valve, and position the balloon within the annular ring. CAUTION: Balloon aortic valvuloplasty of the native aortic valve is required prior to DS insertion to facilitate placement of the THV. CAUTION: For native aortic stenosis, the diameter of the aortic valvuloplasty balloon used for pre-dilatation should be at least equal to the minimum annular diameter. For the ALLEGRA 31, at least a 24 mm balloon is recommended.
7	Pre-dilate the native aortic valve under a short episode of rapid ventricular pacing. If a balloon aortic valvuloplasty is performed within a failing SAV bioprosthesis, the diameter of the balloon should be less than the inner diameter of the failing SAV bioprosthesis. In addition, it is important to ensure that the balloon aortic valvuloplasty catheter is in the central lumen of the failed SAV bioprosthesis. CAUTION: Balloon aortic valvuloplasty should only be performed when the THV is loaded into the DS and ready for implantation. CAUTION: In case of failing SAV bioprosthesis with structural damage, balloon aortic valvuloplasty should be performed with extreme caution.
8	Introduce the flushed DS with the loaded THV over the guidewire under fluoroscopic guidance throughout the introducer sheath to the aorta. Maintain the guidewire position in the left ventricle. NOTE: Before introducing the DS into the patient, visually check the correct loading of the THV.
9	Introduce the DS until the THV crosses the aortic valve. CAUTION: If you feel any resistance during introducing the DS, the procedure must be stopped, and the cause of the resistance has to be evaluated.

9.10. THV DEPLOYMENT

Step	Procedure
1	Position the catheter so that the radiopaque ring of the catheter tip is level with the valve annulus plane marked by the pigtail catheter. In the case of an implantation in a failing aortic valve bioprosthesis, position the catheter to align the THV with the position of the failed SAV bioprosthesis.
2	Hold the DS stable and slowly turn the release handle until the catheter starts to retract and the THV starts to open. CAUTION: Make sure that the catheter at the integrated sheath is stable to avoid any position change.
3	The six radiopaque gold markers of the THV are at the level of the new valve plane and can be used to assess the position in reference to the coronary ostia. NOTE: The upper limit of the THV skirt is also marked by these six radiopaque markers and must not be placed above the coronary ostia.
4	After optimizing the co-axial alignment of the THV to the annular plane and confirming that the new valve plane is below the coronary ostia, continue to rotate the release handle counterclockwise to progressively deploy the THV. Try to maintain the DS in a stable position.
5	Continue to deploy the THV in a controlled manner under fluoroscopic guidance using the radiopaque marker band of the catheter, the radiopaque stent frame, the 6 gold markers and the radiopaque loading anchor to assess valve position. Adjust the position of the THV to achieve a target implant depth of 4 mm below the annular plane. CAUTION: Once the bottom of the THV is implanted in the aortic annulus, do not pull, or push on the delivery system to adjust the implant depth. NOTE: In the bottom of the stent frame, the first half stent cell has a height of 4 mm. NOTE: In case the THV needs to be recaptured, follow the steps of the valve recapturing procedure.
6	If the valve position is acceptable, deploy the THV until the release handle reaches the safety position. CAUTION: After the safety position is reached, do not turn the release handle further. NOTE: The THV is recapturable during deployment before the catheter reaches the safety position. In case the THV must be recaptured, follow the steps of the valve recapturing procedure.
7	Turn the unlock wheel at the safety position by a clockwise movement. NOTE: Make sure that the unlock wheel is engaged after turning. NOTE: Use the lock and unlock symbols on the handle as guidance.
8	Verify an acceptable THV position, pull back the pigtail catheter and fully release the THV by moving the release handle until the THV is fully deployed. CAUTION: After full THV deployment, check if the T-bars are free from the catheter. NOTE: Always pay attention to the stability of the THV.
9	Perform a contrast injection to confirm the position of the THV and the patency of the coronaries. NOTE: All implantation steps should be performed under fluoroscopic control.

9.11. THV RECAPTURE

CAUTION: Recapture and deployment of the THV can be attempted up to three times. If the THV has been fully recaptured three times, it must be removed from the patient and discarded.

NOTE: The THV is recapturable during deployment before the catheter reaches the safety position. Up to this position the T-bars are still attached to the loading anchor.

NOTE: THV recapture is an optional maneuver and should only be undertaken if the initial THV position is not acceptable.

Step	Procedure
1	To recapture the THV rotate the release handle in a clockwise direction. CAUTION: Always perform THV recapture and repositioning procedure only under fluoroscopic guidance. CAUTION: Stop closing the catheter if a resistance occurs. You may check and adjust the angulation between the THV and the cartridge.
2	Depending on the individual case, a partially recaptured THV can be repositioned, or full recapture may be necessary. CAUTION: Always perform THV recapture and repositioning procedure only under fluoroscopic guidance.
3	To fully recapture the THV rotate the release handle until the radiopaque marker band of the cartridge is at the proximal end of the radiopaque tip. CAUTION: Do not over-close the catheter to avoid overlapping. If this occurs the diameter of the tip will be increased by the overlapping catheter and make it difficult to remove the delivery system through the introducer sheath.
4	After recapture (partial or full), carefully reposition the THV to the recommended target depth of 4 mm. CAUTION: A higher implant position may increase the risk for paravalvular leakage. A lower implant position may contribute to a higher risk of conduction disorders.
5	Redeploy the THV by using the release handle.
6	Either complete THV deployment or, if necessary, recapture the THV again. If the THV has been recaptured three times, withdraw the THV from the patient and reassess the situation.

9.12. WITHDRAWAL OF THE DELIVERY SYSTEM

Step	Procedure
1	Retract the guidewire, bringing its soft tip closer to the implanted valve to make this maneuver easier. This will take away the spring force placed by the guidewire against the catheter tip and therefore make the crossing of the THV safer.
2	Slowly withdraw the DS and observe under fluoroscopic guidance if the DS tip is not trapped or dragging the THV, mainly when crossing the T-bars.

Step	Procedure
3	After the tip crosses the implanted valve, gently retract the DS into the descending aorta and close the catheter by turning the release handle clockwise. Close the catheter until the radiopaque marker band of the cartridge is at the proximal end of the radiopaque tip. CAUTION: Do not over-close the catheter to avoid overlapping. If this occurs the diameter of the tip will be increased by the overlapping catheter and make it difficult to remove the delivery system through the introducer sheath. CAUTION: The catheter closure is recommended to be done in a straight portion of the descending aorta. NOTE: The closure allows for withdrawal of the DS out of the introducer sheath without excessive bleeding.

9.13. VERIFICATION OF THV POSITION AND MEASUREMENTS

Step	Procedure
1	Perform a supra-aortic angiogram to evaluate device position, function, and coronary patency. CAUTION: Proceed with caution if it is necessary to post-dilate the THV in the case of unacceptable para-valvular leakage. Post-dilatation could damage the device integrity or cause migration of the THV. NOTE: If there is doubt about the coronary patency, an additional angiogram of each coronary branch can be selectively made using an angiographic catheter.
2	Measure and record the transvalvular pressure gradients.
3	Remove all catheters and sheaths.
4	Perform the closure of the arterial access according to preferred techniques and devices.

10. HOW SUPPLIED

10.1. AVAILABLE THV SIZES

Model ALLEGRA Transcatheter Heart Valve (THV)	ALLEGRA 23	ALLEGRA 27	ALLEGRA 31
Catalogue / Reference number	THV-AO0023G1RE	THV-AO0027G1RE	THV-AO0031G1RE
Frame height [mm]	37.3	41.3	43.0
Inflow diameter [mm]	23.8	27.4	31.0
Outflow diameter [mm]	20.8	24.0	24.0
Sealing height [mm]	12	12	12
Aortic annulus diameter [mm]	16.5 – 21	21 – 24	24 – 27

Aortic annulus perimeter [mm]	51.8 – 66.0	66.0 – 75.4	75.4 – 84.8
Aortic annulus area [mm ²]	214 – 346	346 – 452	452 – 573
Inner diameter aortic valve bioprosthesis [mm]	16.5 – 21	21 – 24	24 – 27

10.2. AVAILABLE DELIVERY SYSTEMS

Model IMPERIA Delivery System	IMPERIA-S	IMPERIA-M
Catalogue / Reference number	DSL-AO18G2115S	DSL-AO18G2115M
Compatible with ALLEGRA THV	ALLEGRA 23	ALLEGRA 27 ALLEGRA 31
Outer diameter	18 Fr (6.0 mm)	18 Fr (6.0 mm)

10.3. STORAGE, PACKAGING AND DISPOSAL OF DEVICES

The THV must be stored between 10°C and 38°C. Avoid locations where extreme temperature fluctuation may occur.

CAUTION: Do not freeze. Always store the THV in a cool, dry place between 10°C and 38°C. Any THV that has been frozen must not be used for human implantation.

The THV is sterilized with and stored in glutaraldehyde solution. The glutaraldehyde storage solution has a concentration of 0.2 %.

CAUTION: Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure or inhalation of the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water. In the event of contact with eyes, seek immediate medical attention.

The THV is supplied sterile in a sealed product jar with a screw cap. The outside of the product jar is non-sterile and must not be placed in the sterile field. A warm and a cold temperature indicator are placed inside the carton box.

CAUTION: Do not use the THV if the warm temperature indicator is black, or the cold temperature indicator's bulb is purple.

NOTE: The product has been exposed to temperatures above 40°C if indicator is black. The product has been exposed to temperatures below 0°C if indicator bulb is purple.

NOTE: After implantation of the THV, the glutaraldehyde solution should be disposed of according to hospital procedure. The package of the THV can be disposed as non-biohazard standard hospital waste.

The DS can be safely stored within temperature limits of 10°C to 30°C and shall be kept dry. The DS and the LS are sterilized with ethylene oxide gas. They are packed together in a double sterile bag. The DS and the LS are sterile if the sterile bags are undamaged and unopened, and the ethylene oxide indicator is green. The outer surface of the outer bag is non-sterile and must not be placed in the sterile field.

CAUTION: Use the DS only if the ethylene oxide indicator on the inner sterile bag is green.

NOTE: The used DS may be disposed of in the same manner that the hospital waste and bio-hazardous materials are handled. There are no special or unusual risks related to the disposal of the DS. The package of the DS can be disposed as standard non-biohazard hospital waste.

11. MRI SAFETY INFORMATION

Nonclinical testing on Magnetic Resonance Imaging (MRI, MR) indicates that the ALLEGRA Transcatheter Heart Valve is MR Conditional. It can be scanned under the following conditions:

- Static magnetic field of ≤ 3 Tesla
- Spatial gradient field of ≤ 1500 Gauss/cm
- Maximum whole-body-averaged Specific Absorption Rate (SAR) of 2.0 W/kg for 15 minutes of scanning.
- Normal mode of MR system

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the THV.

12. INFORMATION ON THE MATERIALS AND SUBSTANCES TO WHICH THE PATIENT CAN BE EXPOSED

Components and materials of the THV to which the patient will be exposed after implantation:

Component	Material	Quantitative composition**
Leaflets, Skirt	Bovine pericardium (cross linked with glutaraldehyde) *	0.96g
Suture material	Polytetrafluoroethylene (PTFE, ePTFE)	0.18g
Support stent structure	Nitinol	1.02g
Radiopaque markers	Gold	0.17g

* As the device contains animal tissue, which is rendered non-viable, it is subject to Commission Regulation (EU) No 722/2012.

** Mass is based on the material of the largest THV size.

Components and materials of the DS:

Component	Material
Catheter System	Polyethylene (PE), Polycarbonate (PC), Polyetheretherketone (PEEK), Polyamide (PA), Polyimide (PI), Polytetrafluoroethylene (PTFE), Polyetherblockamide (PEBA), Polyurethane (PU), Polyoxymethylene (POM), Polyethylene terephthalate (PET), Polyphenylsulfon (PPSU), Thermoplastic polyurethane (TPU)
Hydrophilic coating	Hydrophilic Coating Solution, 10-30% Isopropanol
Radiopaque markers	Platinum/Iridium
Others	Silicone, Stainless steel, Titanium, Rubber, Adhesives

13. PATIENT COUNSELING INFORMATION

Physicians should review the instructions for use for applicable information to be shared with the patient. A patient implant card and a patient leaflet are included in the device package and contains identifying information about the implanted device.

Physicians should provide the patient with information on risks/benefits (e.g., contraindications, warnings, precautions, residual risks, potential undesirable side-effects, and alternative treatments), follow-up visits and medication to be taken. Besides, the physicians should advise the patient to always carry the implant card and direct the patient to the website indicated on the implant card where the patient can find additional information in relation to the implanted device. The implant card contains blank fields to enter the patient's name, date of implantation and healthcare institution filled by handwritten text. Information to identify the implanted device including device name, device model (device type), serial number and Unique Device Identifier (UDI) as well as manufacturer address and website is provided by suitable implant card device stickers. Device stickers are included within the device package of the implantable device. The implant card contains a dedicated field to place the implant card device sticker. The patient leaflet includes information regarding warnings, precautions, or measures, expected lifetime of the device and necessary follow-up information. It also contains information for a healthcare professional with regard to reciprocal interference with reasonably foreseeable external influences, medical examinations or environmental conditions. In addition, the overall qualitative and quantitative information on the materials and substances to which patients can be exposed is included in the patient leaflet.

At some time, the ALLEGRA Transcatheter Heart Valve may need to be replaced. Careful and continuous medical follow-up is advised so that implant related complications can be diagnosed and properly managed. The ALLEGRA Transcatheter Heart Valve was tested to mimic 5 years of use and even tested for 10 years (metal stent frame) according to current applicable standards.

NOTE: The healthcare institution must ensure that the patient receives an implant card completely filled with all the predefined information as well as the information given by the patient leaflet.

NOTE: Only a complete filled implant card allows the patient to have access to device related information.

14. TRAINING

The manufacturer provides detailed loading and operator training. For training material and in-service support, contact a manufacturer representative.

CAUTION: Implantation of the THV may only be performed by physicians who have received the appropriate training.

CAUTION: Loading of the THV may only be performed by physicians or hospital personnel who have received the appropriate training.

15. REPORTING TO MANUFACTURER AND COMPETENT AUTHORITY

Any health care professional (e.g., customer or user of the device) and/or a patient treated with the ALLEGRA THV System who has identified any serious incident that has occurred in relation to the provided device should inform the manufacturer (info-sh@biosensors.com) and the local competent authority immediately. The manufacturer has to be notified by telephone or in written correspondence. When reporting a serious incident provide device name and catalogue / reference number, lot number or serial number and the nature of incident.

16. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE AND EXPECTED CLINICAL BENEFIT

For implantable devices, the manufacturer sets up a regular Summary of Safety and Clinical Performance (SSCP). This report for healthcare professionals and patients will be made publicly available in the future via Eudamed (European Databank on Medical Devices). The SSCP can

be found at <https://ec.europa.eu/tools/eudamed>. Search for the SSCP using the manufacturer and device name, and any of the following elements, as applicable: device model or catalogue / reference number. Until the Eudamed is fully functional the SSCP can be requested at the website of the manufacturer.

The ALLEGRA THV System is designed to relieve aortic stenosis in diseased native aortic valves or to relieve aortic stenosis or aortic insufficiency in failed SAV bioprosthesis. The ALLEGRA Transcatheter Heart Valve is delivered on a catheter based minimal invasive implantation technique, thus eliminating the need for open heart surgery. The ALLEGRA Transcatheter Heart Valve replaces the diseased native or failed SAV bioprosthesis. Successful treatment improves patient activity and quality of life and reduces mortality.

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www.biosensors.com

ALLEGRA and IMPERIA are trademarks or registered trademarks of NVT AG, the holding company of NVT GmbH. ALLEGRA™ and IMPERIA™ are CE approved medical devices. Third party brands are trademarks of their respective intellectual property right owners.

CAUTION: The law restricts these devices to sale by or on the order of a physician and for the use by or under the direction of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device. Not available in the United States and any other country where applicable health authority product registration has not been obtained.

The legal manufacturer of ALLEGRA™ and IMPERIA™ is NVT GmbH, which is part of Blue Sail Medical Co., Ltd. NVT AG and Biosensors International Group, Ltd. and its affiliates are collaborating for the commercialization of the ALLEGRA™ and IMPERIA™ medical devices.

© 2025 NVT GmbH. All rights reserved. Rev.00127502_042025_03_ML

CE
0124

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* Cihaz cansız hâle getirilmiş hayvan dokusu içerdiğinden, 722/2012 Sayılı (AB) Komisyon Yönetmeliğine tabidir.